



Colours Available:



Ultra-Compact Travel Scooter

Enjoy worry-free travel with the Apex Lite, an ultra-compact scooter that is simple to disassemble. Ride in comfort and style with the foldable swivel seat and keep your personal items handy with the front basket.

Specifications	
Scooter	
Weight Capacity	18 st 8 lbs (118 kg)
Overall Length	103 mm (39.5")
Overall Width	495 mm (19.5")
Drive Tyres	Front: 178 mm Solid (2" x 7") Rear: 200 mm Solid (2.5" x 8")
Ground Clearance	Mid-Frame: 570 mm (2.25") At Motor: 340 mm (1.25")
Heaviest Piece When Disassembled	Front Section: 15 kg (33 lbs)
Colour(s)	Red, Blue, Black
Standard Seating	Type Folded Material: Black Vinyl Width: 431.8 mm (17") Depth: 431.8 mm (17") Weight: 1.6 st 23 lbs (10.4 kg)
Overall Height	33.25" - 35.25" (844 - 895 mm)
Mobility	
Turning Radius	112 mm (44.125")
Max Speed	4 mph (6.4 km/h)
Maximum Obstacle Climbing Ability	50 mm (2")
Drive System	Rear - Wheel Drive, Sealed Transaxle, 24v DC Motor
Dual Braking System	Regenerative and Electromechanical
Battery	
Battery Requirements	Standard: 2 x 12 Ah Deep Cycle (Optional) 18 Ah
Battery Charger	2-Amp Off-Board
Battery Weight	Standard: 12 Ah 8.8 kg (19.5 lbs) Optional: 18 Ah 13.4 kg (29.5 lbs)
Total Weight Without Batteries	36.3 kg (80 lbs)

Features



8" Rear Alloy Style Wheels



Easy to Use Scooter Controls



Disassembles into 5 Pieces