



Colours Available:



Ultra-Compact Travel Scooter

Enjoy worry-free travel with the Apex Lite, an ultra-compact scooter that is simple to disassemble. Ride in comfort and style with the foldable swivel seat and keep your personal items handy with the front basket.



Specifications	
Scooter	
Weight Capacity	21 st 6 lbs (136 kg)
Overall Length	1022 mm (40.25")
Overall Width	490 mm (19.25")
Drive Tyres	Front: 203 mm Solid (8") Rear: 203 mm Solid (8")
Ground Clearance	At Motor: 340 mm (1.25")
Heaviest Piece When Disassembled	Front Section: 15 kg (33 lbs)
Colour(s)	Red, Blue, Black
Standard Seating	Type Folded Material: Black Vinyl Width: 431.8 mm (17") Depth: 431.8 mm (17") Weight: 1.6 st 23 lbs (10.4 kg)
Overall Height	33.25" - 35.25" (844 - 895 mm)
Mobility	
Turning Radius	1327 mm (52.25")
Max Speed	4 mph (6.4 km/h)
Maximum Obstacle Climbing Ability	57.5 mm
Drive System	Rear - Wheel Drive, Sealed Transaxle, 24v DC Motor
Dual Braking System	Regenerative and Electromechanical
Battery	
Battery Requirements	2 x 12 Ah
Battery Charger	2-Amp Off-Board
Battery Weight	9 kg
Total Weight Without Batteries	37 kg

Features



8"Rear Alloy Style Wheels



Easy to Use Scooter Controls



Disassembles into 5 Pieces



REV 02 - 30/01/2023

