

MODEL	S19/S19P QUEST
WEIGHT CAPACITY	115kgs(250 lbs.)
SEAT: TYPE/SIZE	16"
DRIVE WHEEL	8"(PU Tire) / 9" (Pneumatic Tire)
FRONT CASTER (WHEEL)	7" (PU Tire) / 8" (Pneumatic Tire)
REAR CASTER (ANTI-TIPPER)	YES
MAX SPEED	6 KM/H
BATTERY SPECIFICATIONS	24V 11.5Ah Li-Polymer / 12V 12AH* 2 pcs
BATTERY RANGE	16 km (Lithium) / 22 Km (Lead acid)
CHARGER TYPE	2.5Amp, 110/240 Volt, Li-ION Charger 2 Amp,120/240 Volt, Lead acid Charger
CONTROLLER TYPE	S-DRIVE 45 Amp
MOTOR TYPE	270W
WEIGHT: W/ BATTERY	26 kgs
WEIGHT : W/O BATTERY	23.1 kgs
Battery (removable, the heaviest part)	2.9 kg(Lithium) /9.2 kg(Lead acid)
TURNING RADIUS	820 mm
LENGTH	970 mm
WIDTH	485 mm
HEIGHT	910~1020 mm
Folded Size	750 mm*500 mm*450mm
SEAT DEPTH	330 mm
BACK HEIGHT	330 mm
WHEEL BASE	695 mm
GROUND CLEARANCE	20 mm
Max Safe Slope	3 Degree
Obstacle Ability	20 mm