

GO CHAIR



Colours Available:



Enhance Your Independence! And Enjoy Your Favourite Places

The Go Chair® is number one for portability! With feather-touch disassembly, you can take your Go Chair to all your favorite places and enjoy greater independence.

Specifications	
Chair	
Weight Capacity	21 st 6 lbs (136 kg)
Overall Length	842 mm (33")
Overall Width	575 mm (22.5")
Drive Tyres	229 mm (9") - Solid
Ground Clearance	40 mm (1.6")
Heaviest Piece When Disassembled	16.5 kg
Colour(s)	Red, Orange, Quartz, Blue, Yellow, White
Power	
Turning Radius	648 mm (33")
Max Speed	Up to 3.8 mph (6 km/h)
Drivetrain	Rear-Wheel Drive, Two Motor
Dual Breaking System	"Intelligent Braking" Electronic Regenerative, Disc Park Brake
Battery	
Battery Requirements	2 x 18 Ah Batteries
Battery Weight	14 06 kg
Battery Charger	18 Ah: 2 - Amp, Off-Board 22 Ah: 3.5 - Amp, Off-Board
Total Weight Without Batteries	44 kg (97 lbs)

Features



Disassembles into Five Pieces



Large Comfortable Foot Plate



Adjustable Height, Width & Armrest

