

GOCHAIR



Colours Available:



Enhance Your Independence! And Enjoy Your Favourite Places

The Go Chair® is number one for portability! With feather-touch disassembly, you can take your Go Chair to all your favorite places and enjoy greater independence.

Features



Disassembles into Five Pieces



Large Comfortable Foot Plate



Adjustable Height, Width & Armrest

Specifications	
Chair	
Weight Capacity	21 st 6 lbs (136 kg)
Overall Length ³	842 mm
Overall Width ³	575 mm
Drive Tyres	228.6 mm (9") - Solid
Ground Clearance ³	40 mm (1.6")
Heaviest Piece When Disassembled ³	16.5 kg
Colour(s)	Ruby Red, Amber Orange, Rose Quartz, Sapphire Blue, Citrine Yellow, Pearl White
Power	
Turning Radius ³	645 mm (25.4")
Max Speed ¹	Up to 3.7 mph
Drivetrain	Rear-Wheel Drive, Two Motor
Dual Breaking System	"Intelligent Braking" Electronic Regenerative, Disc Park Brake
Battery	
Battery Requirements ³	2 x 18 Ah Batteries 2 x 22 Ah Batteries
Battery Weight ⁴	14 kg
Battery Charger	18 Ah: 2 - Amp, Off-Board 22 Ah: 3.5 - Amp, Off-Board
Total Weight Without Batteries ^{3, 5}	30.39 kg

NOTE:

- 1 Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tyre condition. This specification can be subject to a variance of + or -) 10%.
- 2 Includes standard seat weight.
- 3 Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.
- 4 Battery weight may vary based on battery manufacturer.
- 5 Does not meet airline requirements.

All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tyre pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3 %. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.

