

**GO GO**  
**ELITE TRAVELLER®**



Colours Available:



# The Original Travel Scooter

## Make Everyday an Adventure

Every day is a new adventure with the Go-Go Elite Traveller®, an exceptionally stable and maneuverable travel scooter. With easy disassembly and a convenient storage basket, the Go-Go Elite Traveller keeps up with your active life!

# The Original Travel Scooter Make Everyday an Adventure

## Features



Easy Controls



Feather-Touch Disassembly  
With One Hand



Large Back Pocket For  
Extra Storage

Specifications		
Scooter	4-Wheel	3-Wheel
Weight Capacity	21 st 6 lbs (136 kg)	21 st 6 lbs (136 kg)
Overall Length <sup>3</sup>	1003 mm (39.5")	940 mm (37")
Overall Width <sup>3</sup>	495 mm (19.5")	495 mm (19.5")
Drive Tyres	Front: 178 mm (7") - Solid Rear: 203 mm (8") - Solid	Front: 203 mm (8") - Solid Rear: 203 mm (8") - Solid
Ground Clearance <sup>3</sup>	57 mm (1.5")	57 mm (1.5")
Heaviest Piece When Disassembled <sup>3</sup>	16 kg	13.4 kg
Colour(s)	Red & Blue	Red & Blue
<b>Mobility</b>		
Turning Radius <sup>3</sup>	1120 mm (44.1")	838 mm (33")
Max Speed <sup>1</sup>	4 mph (6km/h)	4 mph (6km/h)
<b>Battery</b>		
Battery Requirements <sup>3</sup>	Standard: 2 x 12 Ah Deep Cycle Optional: 2 x 18 Ah Deep Cycle	Standard: 2 x 12 Ah Deep Cycle Optional: 2 x 18 Ah Deep Cycle
Battery Weight <sup>4</sup>	12 Ah - 8.8 kg 18 Ah - 13.4 kg	12 Ah - 8.8 kg 18 Ah - 13.4 kg
Battery Charger	2 - Amp Off-Board	2 - Amp Off-Board
Total Weight Without Batteries <sup>2, 4</sup>	36 kg	34.2 kg

### NOTE:

1 Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tyre condition. This specification can be subject to a variance of (+ or -) 10%.

2 Includes standard seat weight.

3 Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.

4 Battery weight may vary based on battery manufacturer.

All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.

