

GO GO
ELITE TRAVELLER®



Colours Available:



The Original Travel Scooter

Make Everyday an Adventure

Every day is a new adventure with the Go-Go Elite Traveler®, an exceptionally stable and maneuverable travel scooter. With easy disassembly and a convenient storage basket, the Go-Go Elite Traveler keeps up with your active life!

The Original Travel Scooter Make Everyday an Adventure

Specifications		
Scooter	4-Wheel	3-Wheel
Weight Capacity	21 st 6 lbs (136 kg)	21 st 6 lbs (136 kg)
Overall Length	1003 mm (39.5")	940 mm (37")
Overall Width	495 mm (19.5")	495 mm (19.5")
Overall Height	851 - 902 mm (33.5" - 35.5")	851 - 902 mm (33.5" - 35.5")
Drive Tyres	Front: 178 mm (2" x 7") - Solid Rear: 203 mm (2.5" x 8") - Solid	Front: 203 mm (2" x 8") - Solid Rear: 203 mm (2.5" x 8") - Solid
Ground Clearance	35 mm (1.38") at Motor	35 mm (1.38") at Motor
Heaviest Piece When Disassembled	Front Section: 15.5 kg	Front Section: 13.5 kg
Colour(s)	Red & Blue	Red & Blue
Mobility		
Turning Radius	1120 mm (44")	838 mm (33")
Max Speed	4 mph (6.4 km/h)	4 mph (6.4 km/h)
Drive System	Rear-Wheel Drive Sealed Transaxle 24v DC Motor	Rear-Wheel Drive Sealed Transaxle 24v DC Motor
Dual Braking System	Regenerative and Electromechanical	Regenerative and Electromechanical
Maximum Obstacle Climbing Ability	50 mm (2")	50 mm (2")
Battery		
Battery Requirements	2 x 12 Ah Sealed Lead Acid Non Spillable Deep Cycle	2 x 12 Ah Sealed Lead Acid Non Spillable Deep Cycle
Battery Weight	12 Ah - 9 kg Optional: 12 Ah Extended Range Battery Pack 10.75 kg 18 Ah - 13.5 kg Optional: 18 Ah Extended Range Battery Pack 14.5 kg	12 Ah - 9 kg Optional: 12 Ah Extended Range Battery Pack 10.75 kg 18 Ah - 13.5 kg Optional: 18 Ah Extended Range Battery Pack 14.5 kg
Battery Charger	2A Off-Board	2 - Amp Off-Board
Total Weight Without Batteries	36 kg	34.25 kg
Range Per Charge	12 Ah Up to 6 miles (9.75 km) 12 Ah Extended Range Up to 8.2 miles (13.2 km) 18 Ah Up to 9.75 miles (15.5 km) 18 Ah Extended Range Up to 12 miles (19.3 km)	12 Ah Up to 7.25 miles (11.5 km) 12 Ah Extended Range Up to 8.5 miles (13.7 km) 18 Ah Up to 10.5 miles (16.75 km) 18 Ah Extended Range Up to 12.4 miles (20 km)

Features



Easy Controls



Feather-Touch Disassembly With One Hand



Large Back Pocket For Extra Storage