

GO GO
ELITE TRAVELLER®



Colours Available:



The Original Travel Scooter

Make Everyday an Adventure

Every day is a new adventure with the Go Go Elite Traveller®, an exceptionally stable and maneuverable travel scooter. With easy disassembly and a convenient storage basket, the Go Go Elite Traveller keeps up with your active life!

Specifications	
Scooter	
Weight Capacity	21 st 6 lbs (136 kg)
Overall Length	1003 mm (39.5")
Overall Width	495 mm (19.5")
Overall Height	851 - 902 mm (33.5" - 35.5")
Drive Tyres	Front: 178 mm (7") - Solid Rear: 203 mm (8") - Solid
Ground Clearance	35 mm (1.38") at Motor
Heaviest Piece When Disassembled	Front Section: 15. 5 kg
Colour(s)	Red & Blue
Mobility	
Turning Radius	1120 mm (44")
Max Speed	4 mph (6.4 km/h)
Drive System	Rear-Wheel Drive Sealed Transaxle 24v DC Motor
Dual Braking System	Regenerative and Electromechanical
Maximum Obstacle Climbing Ability	50 mm
Battery	
Battery Requirements	2 x 12 Ah (Standard) 2 x 17Ah (Optional)
Battery Weight	12 Ah: 8.75 kg 17 Ah 13.5 kg
Battery Charger	2A Off-Board
Total Weight Without Batteries	36 kg

Features



Easy Controls



Feather-Touch Disassembly With One Hand



Large Back Pocket For Extra Storage



REV 02 - 30/01/2022

