

**GO GO**  
ELITE TRAVELLER LX



Colours Available:



## Sleek Style With Luxurious Features

Every day is a new adventure with the Go-Go Elite Traveler LX<sup>®</sup>, an exceptionally stable and maneuverable travel scooter. With easy disassembly and a convenient storage basket, the Go-Go Elite Traveler LX keeps up with your active life!

## Features



Feather-Touch Disassembly  
With One Hand



Large Back Pocket For Storage



Easy to Use Controls

Specifications	
<b>Scooter</b>	
Weight Capacity	19 st 9 lbs (125 kg)
Overall Length <sup>3</sup>	1003 mm (39.5")
Overall Width <sup>3</sup>	495 mm (19.5")
Drive Tyres	Front: 178 mm (7") - Solid Rear: 203 mm (8") - Solid
Ground Clearance <sup>3</sup>	63.5 mm (1.5")
Heaviest Piece When Disassembled <sup>3</sup>	15.2 kg
Suspension	Comfort-Trac Suspension
Colour(s)	Red & blue
<b>Mobility</b>	
Turning Radius <sup>3</sup>	1385 mm (54.5")
Max Speed <sup>1</sup>	4 mph (6 km/h)
<b>Battery</b>	
Battery Requirements <sup>3</sup>	2 x 12 Ah Deep Cycle (Standard) 2 x 18 Ah Deep Cycle (Optional)
Battery Weight <sup>4</sup>	12 Ah - 8.8 kg 18 Ah - 13.5 kg
Battery Charger	2 - Amp Off-Board
Total Weight Without Batteries <sup>3, 4</sup>	35.8 kg

### NOTE:

- 1 Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tyre condition. This specification can be subject to a variance of + or - 10%.
- 2 Includes standard seat weight.
- 3 Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.
- 4 Battery weight may vary based on battery manufacturer.
- 5 Does not meet airline requirements.

All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.

