

GO GO
ELITE TRAVELLER® PLUS



Colours Available:



The Quickest Way To Travel

Perfect For Travelling

It's perfect for travelling, PLUS so much more! The Go-Go Elite Traveller® Plus offers easy disassembly and includes an ergonomic tiller for one-handed operation. Today's fast-paced world is no match for this stylish and portable scooter.

The Quickest Way To Travel Perfect For Travelling

Features



Easy to Use and Straight Forward Controls



Available in both 3 & 4 Wheel

Specifications		
Scooter	4-Wheel	3-Wheel
Weight Capacity	23 st 3 lbs (147 kg)	23 st 3 lbs (147 kg)
Overall Length ³	1023 mm (40.25")	1016 mm (40")
Overall Width ³	540 mm (21.25")	540 mm (21.25")
Drive Tyres	Front: 203 mm (8") - Solid Rear: 229 mm (9") - Solid	Front: 203 mm (8") - Solid Rear: 229 mm (9") - Solid
Ground Clearance ³	38 mm (1.5")	38 mm (1.5")
Heaviest Piece When Disassembled ³	16 kg	15 kg
Colour(s)	Red & Blue	Red & Blue
Mobility		
Turning Radius ³	1162 mm (45.75")	876 mm (34.5")
Max Speed ¹	4 mph (6 km/h)	4 mph (6 km/h)
Battery		
Battery Requirements ³	2 x 18 Ah Deep Cycle	2 x 18 Ah Deep Cycle
Battery Weight ⁴	18 Ah - 13.4 kg	18 Ah - 13.4 kg
Battery Charger	2 - Amp Off-Board	2 - Amp Off-Board
Total Weight Without Batteries ^{3, 4}	44 kg	39.5 kg

NOTE:

1 Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tyre condition. This specification can be subject to a variance of (+ or -) 10%.

2 Includes standard seat weight.

3 Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.

4 Battery weight may vary based on battery manufacturer.

All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3 %. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.



PRIDE
MOBILITY PRODUCTS LTD.
LIVE YOUR BEST®

www.pride-mobility.co.uk