

Jazzy air₂



Colours Available:



Elevate Your Life To The Next Level!

At the touch of a switch, the Jazzy Air[®] 2 elevates 12 inches in just 11 seconds. With a sleek, modern look, this power chair safely drives at speeds up to 4 mph while elevated, allowing users to enjoy face-to-face social engagement.

Specifications	
Chair	
Weight Capacity	21 st 6 lbs (136 kg)
Overall Length	1067 mm (42")
Overall Width	590 mm (23.25")
Front Tyres Rear Tyres	Front: 152mm (6") Solid Rear: 254 mm (10") Solid
Ground Clearance	25 mm (1")
Suspension	Active-Trac Suspension With Mid-Wheel 6 Technology
Colour(s)	Citrine yellow, Pearly white, Ruby red, Sapphire blue, Black pearl, Garnet red, Black onyx, Pink topaz, Silver tanzanite
Heaviest Piece When Disassembled	82 kg
Power	
Turning Radius	616 mm (24.25")
Max Speed	4 mph (6.4 km/h)
Battery	
Battery Requirements	2 x 50 Ah, Deep Cycle
Battery Weight	14.5 kg
Battery Charger	5 - Amp Off-Board
Total Weight Without Batteries	91.5 kg


Features



Fold Flat Seat



Drives Safely at 4mph While Elevated at 12"



Angle, and Height-Adjustable Headrest



REV 02 - 30/01/2023

