

Jazzy air₂



Colours Available:



Elevate Your Life To The Next Level!

At the touch of a switch, the Jazzy Air[®] 2 elevates 12 inches in just 11 seconds. With a sleek, modern look, this power chair safely drives at speeds up to 4 mph while elevated, allowing users to enjoy face-to-face social engagement.

| Specifications | |
|----------------------------------|--|
| Chair | |
| Weight Capacity | 21 st 6 lbs (136 kg) |
| Overall Length | 1067 mm (41,93") |
| Overall Width | 590 mm (23.19") |
| Drive Tyres | Drive: 254 mm (10") Solid Tyres |
| Ground Clearance | 25 mm (1") |
| Suspension | Active-Trac Suspension With Mid-Wheel 6 Technology |
| Colour(s) | Citrine yellow, Pearly white, Ruby red, Sapphire blue, Black pearl, Garnet red, Black onyx, Pink topaz, Silver tanzanite |
| Heaviest Piece When Disassembled | 82 kg |
| Power | |
| Turning Radius | 616 mm (24.25") |
| Max Speed | 4 mph (6 km/h) |
| Battery | |
| Battery Requirements | 2 x 50 Ah, Deep Cycle |
| Battery Weight | 16.9 kg |
| Battery Charger | 5 - Amp Off-Board |
| Total Weight Without Batteries | 81.5 kg (180 lbs) |


Features



Fold Flat Seat



Drives Safely at 4mph While Elevated at 12"



Angle, and Height-Adjustable Headrest