

QUEST



Colours Available:



Light, Portable And Easy To Fold

With the Pride Quest Folding Scooter, the world is at your fingertips! With a unique folding design, this scooter folds in just three simple steps and fits easily in a car boot. Enjoy full lighting and speeds up to 4 mph.

| Specifications | |
|--------------------------------|--|
| Scooter | |
| Weight Capacity | 17 st 10 lbs (115 kg) |
| Overall Length | 970 mm (38") |
| Overall Width | 485 mm (19") |
| Overall Height | 910 - 1020 mm (35.8" - 40.15) |
| Drive Tyres | 203 mm (8") PU Tyre 228 mm (9") Pneumatic Tyre |
| Ground Clearance | 20 mm |
| Colour(s) | Black, Red |
| Suspension | Indego Suspension |
| Folded Size | 750 mm 500 mm 450 mm |
| Mobility | |
| Turning Radius | 820 mm (32.2") |
| Max Speed | 4 mph (6.4km/h) |
| Maximum Obstacle Climbing | 20 mm |
| Drivetrain | 4-Pole 270w |
| Battery | |
| Battery Requirements | 11.5 Ah Li |
| Battery Weight | 3 kg |
| Battery Charger | 2.5 - Amp 110/240v, Li- ION Charger 2 Amp, 120/240v Lead Acid Charger |
| Total Weight Without Batteries | 23 kg |

Features



Easy Controls



Front and Rear LED Lighting



Easy to Transport



REV 02 - 30/01/2023

