


# QUEST



Colours Available:



## Light, Portable And Easy To Fold

With the Pride Quest Folding Scooter, the world is at your fingertips! With a unique folding design, this scooter folds in just three simple steps and fits easily in a car boot. Enjoy full lighting and speeds up to 4 mph.

## Features



Easy Controls



Front and Rear LED Lighting



Easy to Transport

Specifications	
<b>Scooter</b>	
Weight Capacity	18 st (115 kg)
Overall Length <sup>3</sup>	930 mm (36.6")
Overall Width <sup>3</sup>	485 mm (19.09")
Drive Tyres	203mm
Ground Clearance <sup>3</sup>	60 mm
Colour(s)	Red, Black, White
Suspension	Indego Suspension
<b>Mobility</b>	
Turning Radius <sup>3</sup>	820 mm (32.2")
Max Speed <sup>1</sup>	4 mph (6 km/h)
Max. Climb Angle	7 Degrees
Drivetrain	4-Pole 270w
<b>Battery</b>	
Battery Requirements <sup>3</sup>	2 x 11.5 Ah Battery
Battery Weight <sup>4</sup>	3.6 kg
Battery Charger	2.5 - Amp Off-Board
Total Weight Without Batteries <sup>3,4</sup>	22 kg

### NOTE:

- 1 Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tyre condition. This specification can be subject to a variance of + or -) 10%.
- 2 Includes standard seat weight.
- 3 Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.
- 4 Battery weight may vary based on battery manufacturer.
- 5 Does not meet airline requirements.

All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.

