

REVO 2.0



Colours Available:



# Robustly Rugged And Perfectly Convenient

Robust meets convenience. With the Revo® 2.0, there is nothing you can't handle! This tough and dependable scooter has powerful suspension for a smooth ride, yet easily disassembles for full portability.

# Robustly Rugged And Perfectly Convenient

Specifications	
<b>Scooter</b>	
Weight Capacity	28 st 8 lbs (181 kg)
Overall Length	1194 mm (47")
Overall Width	525 mm (20.75")
Overall Height	955 mm (38")
Drive Tyres	Front: 229 mm (9") Solid, Pneumatic Rear: 229 mm (9") Solid, Pneumatic
Ground Clearance	32 mm (1.25") at Motor
Colour(s)	True Blue, Grey Street
Suspension	Comfort Trac
Heaviest Piece When Disassembled	25.85 kg
<b>Mobility</b>	
Drive System	Rear-Wheel Drive, Sealed Transaxle, 24v DC Motor
Dual Braking System	Regenerative and Electromechanical
Turning Radius	1323 mm (52")
Max Speed	4 mph (6.4km/h)
Maximum Obstacle Climb	65 mm
<b>Battery</b>	
Battery Requirements	2 x 30-35 Ah
Battery Weight	U1 Battery Pack 11.25 kg Each
Battery Charger	3.5A Off-Board
Total Weight Without Batteries	60 kg

## Features



Easily Disassembled Into Five Pieces



Underseat Storage



Comfort Trac Front and Rear Suspension



Full LED Lights Package



REV 02 - 30/01/2023

**PRIDE**  
MOBILITY PRODUCTS LTD.  
LIVE YOUR BEST®

[www.pride-mobility.co.uk](http://www.pride-mobility.co.uk)