

**REVO 2.0**



Colours Available:



# Robustly Rugged And Perfectly Convenient

Robust meets convenience. With the Revo® 2.0, there is nothing you can't handle! This tough and dependable scooter has powerful suspension for a smooth ride, yet easily disassembles for full portability.

# Robustly Rugged And Perfectly Convenient

## Features



Easily Disassembled Into Five Pieces



Underseat Storage



Comfort Trac Front and Rear Suspension



Full LED Lights Package

| Specifications                                 |                        |
|--|------------------------|
| <b>Scooter</b>                                 |                        |
| Weight Capacity                                | 26 st 11 lbs (170 kg)  |
| Overall Length <sup>3</sup>                    | 1194 mm (47")          |
| Overall Width <sup>3</sup>                     | 525 mm                 |
| Drive Tyres                                    | 228.5 mm               |
| Ground Clearance <sup>3</sup>                  | 32 mm (1.25") at Motor |
| Colour(s)                                      | Grey Street, True Blue |
| Suspension                                     | Comfort Trac           |
| Heaviest Piece When Disassembled <sup>3</sup>  | 25.85 kg               |
| <b>Mobility</b>                                |                        |
| Turning Radius <sup>3</sup>                    | 1321 mm (52")          |
| Max Speed <sup>1</sup>                         | 4 mph (6 km/h)         |
| <b>Battery</b>                                 |                        |
| Battery Requirements <sup>3</sup>              | 2 x 30 Ah Batteries    |
| Battery Weight <sup>4</sup>                    | 11.5 kg                |
| Battery Charger                                | 3.5 - Amp Off-Board    |
| Total Weight Without Batteries <sup>3, 4</sup> | 60 kg                  |

### NOTE:

- 1 Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tyre condition. This specification can be subject to a variance of + or - 10%.
- 2 Includes standard seat weight.
- 3 Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.
- 4 Battery weight may vary based on battery manufacturer.
- 5 Does not meet airline requirements.

All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3 %. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.



**PRIDE**  
MOBILITY PRODUCTS LTD.  
LIVE YOUR BEST®

[www.pride-mobility.co.uk](http://www.pride-mobility.co.uk)