

**ZERO
TURN⁸**

iTURN[™]
TECHNOLOGY



Colours Available:



Exceptional Turning Power With a Incredible Design

Your mobility has turned a corner with the Zero Turn 8, an innovative scooter equipped with iTurn Technology[™], enabling you to maneuver this scooter through small spaces and around tight corners! Plus, enjoy full lighting, plenty of storage and a USB charger.

Exceptional Turning Power With a Incredible Design



Features



Front and Rear LED Lighting



Exclusive & Stylish Low Profile Tyres



Comfort-Trac Suspension



Convenient Front Basket

Specifications	
Scooter	
Weight Capacity	23 st 3 lbs (147 kg)
Overall Length ³	1162 mm (45.75")
Overall Width ³	525 mm (20.75")
Drive Tyres	Front: 178 mm (7") - Solid Rear: 229 mm (9") - Solid
Ground Clearance ³	33 mm (1.3") at Motor
Colour(s)	Fire Opal, Onyx Black, Pearl White, Pewter, Sapphire Blue, Green Moonstone
Suspension	Comfort Trac
Heaviest Piece When Disassembled ³	21.7 kg
Mobility	
Turning Radius ³	948 mm (37.25")
Max Speed ¹	4 mph (6 km/h)
Battery	
Battery Requirements ³	2 x 18 Ah Batteries
Battery Weight ⁴	15 kg
Battery Charger	3.5 - Amp Off-Board
Total Weight Without Batteries ^{3, 4}	54 kg

NOTE:

- 1 Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tyre condition. This specification can be subject to a variance of + or - 10%.
- 2 Includes standard seat weight.
- 3 Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.
- 4 Battery weight may vary based on battery manufacturer.
- 5 Does not meet airline requirements.

All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3 %. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.



PRIDE
MOBILITY PRODUCTS LTD.
LIVE YOUR BEST®