

**ZERO
TURN⁸**

iTURN[™]
TECHNOLOGY



Colours Available:




Exceptional Turning Power With a Incredible Design

Your mobility has turned a corner with the Zero Turn 8, an innovative scooter equipped with iTurn Technology™, enabling you to maneuver this scooter through small spaces and around tight corners! Plus, enjoy full lighting, plenty of storage and a USB charger.

Specifications	
Scooter	
Weight Capacity	23 st 3 lbs (147 kg)
Overall Length	1162 mm (45.75")
Overall Width	525 mm (20.75")
Overall Height	883 mm - 925 mm (34.75" - 36.25")
Drive Tyres	Front: 178 mm (7") - Solid Rear: 229 mm (9") - Solid
Ground Clearance	33 mm (1.3") at Motor
Colour(s)	Fire Opal, Onyx Black, Pearl White, Pewter, Sapphire Blue, Green Moonstone
Suspension	Comfort Trac
Mobility	
Drive System	Rear - Wheel Drive, Two Motor, 24v DC Motor
Dual Braking System	Regenerative and Electromechanical
Turning Radius	948 mm (37.25")
Max Speed	4 mph
Maximum Obstacle Climbing Ability	No Run-Up: 30 mm (1.25") 500 mm Run-Up 50 mm (2")
Battery	
Battery Requirements	2 x 22 Ah 12v Sealed Lead-Acid, Non - Spillable, Deep-Cycle
Battery Weight	22 Ah Battery Pack 15 kg (32 lbs)
Battery Charger	3.5 A Off-Board
Total Weight Without Batteries	54 kg

Features



Front and Rear LED Lighting



Exclusive & Stylish Low Profile Tyres



Comfort-Trac Suspension



Convenient Front Basket