



Colours Available:



The Ultimate Travel Chair

It's Time to Make Your Mark!

The Go-Go® Elite Traveller Sport offers simple disassembly and an ergonomic delta tiller, perfect for one-handed operation. Enjoy the versatility of this sporty scooter.

Features



Compact
Lightweight Design



Simple One Touch
Folding System



Micron Control System

Specifications

Specifications	
Chair	
Weight Capacity	18 st 8 lb (120 kg)
Overall Length ³	876 mm
Overall Width ³	600 mm (23.6")
Drive Tyres	210.6 mm
Ground Clearance ³	44 mm
Power	
Turning Radius ³	900 mm (35.4")
Max Speed ¹	3.7 mph (5.2 km/h)
Battery	
Battery Requirements ³	2 x 6.6 Ah Lithium Battery
Battery Weight ⁴	1.35 kg
Battery Charger	5 - Amp, Off-Board
Total Weight Without Batteries ^{2, 4}	19.8 kg

NOTE:

- 1 Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tyre condition. This specification can be subject to a variance of + or -) 10%.
- 2 Includes standard seat weight.
- 3 Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.
- 4 Battery weight may vary based on battery manufacturer.
- 5 Does not meet airline requirements.

All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tyre pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3 %. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.