



Colours Available:



The Ultimate Travel Chair It's Time to Make Your Mark!

This compact powerchair has plenty of storage and folds in just a few steps to make airport check-in effortless. When your plane touches down at your destination, you are ready to experience the adventures that await you!

Features



Compact
Lightweight Design



Simple One Touch
Folding System



Micron Control System

Specifications	
Chair	
Weight Capacity	18 st 12 lb (120 kg)
Overall Length ³	900 mm (35.4")
Overall Width ³	600 mm (23.6")
Drive Tyres	Front 178mm (7") Rear 216 mm (8.5")
Ground Clearance ³	44 mm
Colour(s)	Silver
Power	
Turning Radius ³	900 mm (35.4")
Max Speed ¹	3 mph (4.83 km/h)
Battery	
Battery Requirements ³	2 x 6.6 Ah Lithium (standard) 2 x 10 Ah Lithium (optional)
Battery Weight ⁴	6.6 Ah: 1.4 kg
Battery Charger	5 - Amp, Off-Board
Total Weight Without Batteries ^{2, 4}	6.6 Ah: 1.4 kg 10 Ah: 2.5 kg

NOTE:

- 1 Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tyre condition. This specification can be subject to a variance of + or - 10%.
- 2 Includes standard seat weight.
- 3 Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.
- 4 Battery weight may vary based on battery manufacturer.
- 5 Does not meet airline requirements.

All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tyre pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.