

**ZERO
TURN 10**

iTURN
TECHNOLOGY



Colours Available:



Regain Your Independence!

When Your Imagination Becomes Reality

The ZT10 combines powerful dual motors and two-wheel drive with revolutionary iturn technology. All with a weight capacity of 28st 8lbs and speeds of up to 6 mph.

Specifications	
Scooter	
Weight Capacity	28 st 8 lbs (182 kg)
Overall Length	1219 mm (48")
Overall Width	613 mm (24.25")
Overall Height	1200 mm - 1238 mm (47.25" - 48.75")
Drive Tyres	Front: 229 mm (9") - Solid Rear: 273 mm (10.75") - Solid
Ground Clearance	51 mm (2") at Motor
Colour(s)	White, Red, Blue
Suspension	Comfort Trac
Heaviest Piece When Disassembled	34 kg
Mobility	
Drive System	Rear - Wheel Drive, Dual Inline, 24v DC Motor
Dual Braking System	Regenerative and Electromechanical
Turning Radius	1092 mm (43")
Max Speed	6 mph (9.6 km/h)
Maximum Obstacle Climbing Ability	75 mm
Battery	
Battery Requirements	2 x 40 Ah
Battery Weight	13.75 kg
Battery Charger	3.5A Off-Board
Total Weight Without Batteries	82.5 kg

Features



Exquisite Front and Rear Light Package



Exclusive iTurn Technology



10.5" Alloy Style Tyres