

iGO FOLD



Colours Available:



Remote Controlled Automatic Folding Power Chair

Experience the convenience of a remote controlled automatic folding powerchair. This compact powerchair folds in just a few seconds with a simple press of a button for effortless transport and storage. With convenient storage pockets, the iGo Fold makes travelling a breeze.

Specifications	
Chair	
Weight Capacity	18 st 8 lb (120 kg)
Overall Length	1100 mm (43.3")
Overall Width	620 mm (24.4")
Drive Tyres	Front: 195 mm (7.6") Rear: 250 mm (9.8")
Ground Clearance	40 mm
Suspension	Limited
Colour(s)	Silver & Black
Power	
Turning Radius	650 mm (25.59")
Max Speed	5 mph (8 km/h)
Drivetrain	2 x 24v / 200w Brushed Motors
Dual Breaking System	Electromechanical
Battery	
Battery Requirements	1 x 20 Ah Lithium Battery
Battery Weight	4.6 kg
Battery Charger	5 - Amp, Off-Board
Total Weight Without Batteries	27.8 kg

Features



Compact
Lightweight Design



Remote Controlled Automatic
Folding System



Back Seat Pocket For
Extra Storage While on-the-go



REV 03 - 09/02/2023

