

**iGO FOLD**



Colours Available:



# Remote Controlled Automatic Folding Power Chair

Experience the convenience of a remote controlled automatic folding power chair. This compact power chair folds in just a few seconds with a simple press of a button for effortless transport and storage. With convenient storage pockets, the iGo Fold makes traveling a breeze.

## Features



Compact  
Lightweight Design



Remote Controlled Automatic  
Folding System



Back Seat Pocket For  
Extra Storage While on-the-go

Specifications	
<b>Chair</b>	
Weight Capacity	18 st 8 lb (120 kg)
Overall Length <sup>3</sup>	1100 mm (43.3")
Overall Width <sup>3</sup>	620 mm (24.4")
Drive Tyres	250 x 60 mm Solid
Ground Clearance <sup>3</sup>	40 mm at Motor
Castors	195 x 45 mm Solid
Suspension	Limited
Colour(s)	Silver & Black
<b>Power</b>	
Turning Radius <sup>3</sup>	650 mm (25.59")
Max Speed <sup>1</sup>	Up to 4 mph
Drivetrain	2 x 24v / 200w Brushed Motors
Dual Breaking System	Electromechanical
<b>Battery</b>	
Battery Requirements <sup>3</sup>	12 Ah Lithium Battery
Battery Weight <sup>4</sup>	4.6 kg
Battery Charger	5 - Amp, Off-Board
Total Weight Without Batteries <sup>3, 4</sup>	27.8 kg

### NOTE:

- 1 Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tyre condition. This specification can be subject to a variance of + or -) 10%.
- 2 Includes standard seat weight.
- 3 Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.
- 4 Battery weight may vary based on battery manufacturer.
- 5 Does not meet airline requirements.

All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3 %. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.

