

iGO lite



Colours Available:



Carbon Fibre, Ultra Lightweight Simply Unfold And Go

Sleek, practical and ultra-lightweight. The all new iGo Lite weighs only 19.8 kg making it our lightest power chair yet! With a built-in USB charger, you'll be able to travel and stay connected with the world around you. It's flight-safe batteries mean the sky is your limit - simply unfold and go, for an adventure of a lifetime.

| Specifications | |
|--------------------------------|---|
| Chair | |
| Weight Capacity | 21 st 3 lbs (136 kg) |
| Overall Length | 960 mm (37.79") |
| Overall Width | 610 mm (24") |
| Overall Weight | 19.8 kg (43.6 lbs) |
| Drive Tyres | Front: 178 mm (7") Rear: 216 mm (8.5") |
| Castors | 195 x 45 mm Solid |
| Suspension | Limited |
| Colour(s) | Black |
| Front Wheel | 178 mm (7") - Solid |
| Power | |
| Turning Radius | 940 mm (35.43") |
| Max Speed | Up to 3.7 mph (6 km/h) |
| Battery | |
| Battery Requirements | 1 x 12 Ah Lithium-Ion Battery |
| Battery Weight | 1.8 kg |
| Battery Charger | 2 - Amp, On Board and Off-Board |
| Total Weight Without Batteries | 18 kg (39.6 lbs) |

Features



Compact
Lightweight Design



On-Board USB Charger



Large Secure Under seat
Storage Compartment While
on-the-go