

iGO lite



Colours Available:



Carbon Fibre, Ultra Lightweight Simply Unfold And Go

The iGo Lite weighs only 18 kg making it our lightest power chair yet! With a built-in USB charger, you'll be able to travel and stay connected with the world around you. It's flight-safe batteries mean the sky is your limit - simply unfold and go, for an adventure of a lifetime.

Specifications	
Chair	
Weight Capacity	21 st 3 lbs (135 kg)
Overall Length	940 mm (37")
Overall Width	610 mm (24")
Drive Tyres	Front: 178 mm (7") solid Rear: 216 mm (8.5") solid
Colour(s)	Black carbon fibre
Power	
Turning Radius	900 mm (35")
Max Speed	4 mph (6.4 km/h)
Battery	
Battery Requirements	1 x 12 Ah Lithium-Ion Battery
Battery Weight	1.8 kg
Battery Charger	2 - Amp, On Board and Off-Board
Total Weight Without Batteries	18 kg

Features



Compact Lightweight Design



On-Board USB Charger



Large Secure Under seat Storage Compartment While on-the-go

UK CA | CE

REV 02 - 30/01/2023

PRIDE[®]
MOBILITY PRODUCTS LTD.
LIVE YOUR BEST[®]