

iGO lite



Colours Available:



Carbon Fibre, Ultra Lightweight Simply Unfold And Go

Sleek, practical and ultra-lightweight. The all new iGo Lite weighs only 19.8 kg making it our lightest power chair yet! With a built-in USB charger, you'll be able to travel and stay connected with the world around you. It's flight-safe batteries mean the sky is your limit - simply unfold and go, for an adventure of a lifetime.

Features



Compact
Lightweight Design



On-Board USB Charger



Large Secure Under seat
Storage Compartment While
on-the-go

Specifications

Specifications	
Chair	
Weight Capacity	21 st 3 lbs (135 kg)
Overall Length ³	960 mm (37.79")
Overall Width ³	610 mm (24")
Overall Weight	19.8 kg (43.6 lbs)
Drive Tyres	216 mm (8.5") - Solid
Castors	195 x 45 mm Solid
Suspension	Limited
Colour(s)	Black
Front Wheel	178 mm (7") - Solid
Power	
Turning Radius ³	900 mm (35.43")
Max Speed ¹	Up to 3.7 mph
Battery	
Battery Requirements ³	12 Ah Lithium-Ion Battery
Battery Weight ⁴	1.8 kg
Battery Charger	2 - Amp, On Board and Off-Board
Total Weight Without Batteries ^{2, 4}	18 kg

NOTE:

- 1 Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tyre condition. This specification can be subject to a variance of (+ or -) 10%.
- 2 Includes standard seat weight.
- 3 Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.
- 4 Battery weight may vary based on battery manufacturer.

All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.