

iRIDE™



Colours Available:



Lightweight Design, For Effortless Travel


Journey far and wide with the iRide! This ultra-lightweight scooter is super portable and perfect for traveling. Equipped with both a manual hand brake and supplementary foot brake, the iRide delivers a safe and comfortable ride!

Specifications	
Scooter	
Weight Capacity	17 st 12 lbs (113 kg)
Overall Length	968 mm (38")
Overall Width	485 mm (19.25")
Overall Height	851 mm (33.5")
Drive Tyres	Front: 179 mm (7") - Solid Rear: 168 mm (6.5") - Solid
Ground Clearance	33 mm
Colour(s)	Orange, Black, Red, Green, Blue
Standard Seating	Type: Compact Material: Molded Plastic Width: 375 mm (14.75") Depth: 343 mm (13.5") Weight: 4.5 kg (10 lbs)
Mobility	
Turning Radius	851 mm (33.5")
Max Speed	Up to 4mph (6.4 km/h)
Drivetrain	36 v Brushless Hub Motor
Dual Breaking System	Dual Disc Hand Brake, Supplemental Foot Brake
Maximum Obstacle Climbing Ability	40 mm
Battery	
Battery Requirements	5.2 Ah
Battery Weight	2 kg
Battery Charger	36v Lithium Off-Board Charger
Total Weight Without Batteries	19 kg

Features



Fordable Tiller With Adjustment Knob



Simple and Easy to Use Controls



Portable and Easy to Travel



REV 02 - 30/01/2023

