

iRIDE™



Colours Available:



Lightweight Design, For Effortless Travel

Journey far and wide with the iRide! This ultra-lightweight scooter is super portable and perfect for travelling. Equipped with both a manual hand brake and supplementary foot brake, the iRide delivers a safe and comfortable ride!


Features



Fordable Tiller With Adjustment Knob



Portable and Easy to Travel



Simple and Easy to Use Controls

Specifications	
Scooter	
Weight Capacity	17 st 12 lbs (113 kg)
Overall Length ³	968 mm (38")
Overall Width ³	488 mm (19.25")
Drive Tyres	Front: 179 mm (7") - Solid Rear: 168 mm (6.5") - Solid
Ground Clearance ³	33 mm (1.3")
Colour(s)	Lime, Mango, Blackberry, Arctic Ice, Raspberry
Seating	Compact
Mobility	
Turning Radius ³	851 mm (33.5")
Max Speed ¹	Up to 3.7 mph (6 km / h)
Drivetrain	36v Brushless Hub Motor
Dual Breaking System	Dual Disc Hand Brake, Supplementary Foot Brake
Battery	
Battery Requirements ³	5.2 Ah Lithium Battery
Battery Weight ⁴	2 kg
Battery Charger	36v Lithium Off-Board Charger
Total Weight Without Batteries ^{2,3}	19 kg

NOTE:

- 1 Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tyre condition. This specification can be subject to a variance of + or -) 10%.
- 2 Includes standard seat weight.
- 3 Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.
- 4 Battery weight may vary based on battery manufacturer.
- 5 Does not meet airline requirements.

All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3 %. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.

