



 The
Hudson
Range

Viva LIFT!™
POWER RECLINERS
By Pride® Mobility Products Corp.





The exquisite Hudson Range from Pride Mobility is proudly manufactured and expertly hand crafted within the UK. The range offers sophisticated design, expert craftsmanship and a variety of comfort and style options that can fit in any home.

Recognised as practical and purpose built chairs with comfort and support in mind, that are engineered for outstanding construction and long-lasting durability. The Hudson Range includes unique features designed to increase comfort, support, quality and reliability.

We are happy to offer a 5 Year Warranty on Electronics and Mechanisms.

PRIDE[®]
MOBILITY PRODUCTS LTD.
LIVE YOUR BEST[®]

CHOOSE FROM ONE OF OUR FOUR MECHANISMS



EXTENDING FOOTREST



POWERED HEADREST

UK
CA | C €

UPWARD LIFT

PREMIER RISE RECLINER CHAIR

The power adjustable headrest adjusts the position of your head, neck and shoulders for maximum support. The power adjustable lumbar enables you to personalize your comfort. The footrest extension allows you to relax and stretch out completely. The USB remote adjusts all functions independently and allows you to charge your smart device. Complete with a safety power cut off switch and emergency single use battery backup. Left or right hand side pockets. For ease of movement these chairs are available with castors or glides. This mechanism has 4 individual motors.

DELUXE RISE RECLINER CHAIR

The footrest extension allows you to relax and stretch out completely. The USB remote adjusts all functions independently and allows you to charge your smart device. Complete with a safety power cut off switch and emergency single use battery backup. Left or right hand side pockets. For ease of movement these chairs are available with castors or glides. This mechanism has 2 individual motors.

TILT IN SPACE DUAL MOTOR

The smooth dual operation chairs quite simply offer independence of the back from the footrest, giving you the option to find that perfect sitting or relaxing position. In addition to this, it can drop your hips below your feet and ankles which is excellent for circulation and water retention.

TILT IN SPACE SINGLE MOTOR

The very best in British designed technology, this mechanism provides a wonderfully smooth recline and is not like any other "Single Motor" chairs. The back angle to seat never changes (fixed) on the recline providing no shear on the back and is ideal for anyone who suffers with back & hip problems. It also comes with a higher foot elevation, again good for circulation.



Optional Chair Backs



Letterbox



Lateral




Waterfall


Standard Fabric Options



Vancouver 21
French Grey fr



Vancouver 21
Mocha fr



Vancouver 21
Praline fr



BUXTON

The Buxton is a great choice of chair if you are seeking something a little more modern looking, whilst still bringing you superior relaxation and comfort. Our exclusive and uniquely designed arms are shaped to add a clean look and finish to this classic. You'll always find your most comfortable and supportive position. It's available, as all of our chairs, with all 4 mechanisms, backs and accessory choices.





DORCHESTER



Optional Chair Backs



Letterbox




Lateral




Waterfall


Standard Fabric Options



Vancouver 21
French Grey fr



Vancouver 21
Mocha fr



Vancouver 21
Praline fr

The ultimate in riser recliner chairs. The Dorchester is our award-winning armchair. Superbly crafted upholstery is complemented by its traditional, classic and stylish looks, with scrolled arms and elegant wings - The Dorchester is everything you could want in a riser recliner. Find your most comfortable position at the touch of the button.

Available with the popular waterfall back, lateral support cushion back and letterbox back. This model is available in all 4 of our mechanisms to suit all of your needs.



The Hereford chair exudes class and style with our exclusive wooden knuckle arms which are complemented by a sweeping curved outer arm. The Hereford is a great choice if you struggle to get in and out of your chair as the knuckles provide you with a good grip whilst the chair is rising. The perfect addition to your living room or office. A truly elegant, sophisticated recliner.

Available in all 4 of our mechanisms, all 3 backs and with 3 options of quality wood finishes: Natural, Teak and Mahogany.

HEREFORD



Optional Chair Backs



Letterbox




Lateral




Waterfall


Standard Fabric Options



Vancouver 21
French Grey fr



Vancouver 21
Mocha fr



Vancouver 21
Praline fr



MONMOUTH

Optional Chair Backs



Letterbox



Lateral




Waterfall


Standard Fabric Options



Vancouver 21
French Grey fr



Vancouver 21
Mocha fr



Vancouver 21
Praline fr



Designed with elegant and subtle scroll arms, the Monmouth is similar to the Hereford in style but without the wooden knuckles, giving you a modern but still classic look. Its classic sweep outer arm gives the chair an aesthetically pleasing look and delivers comfort and relaxation. A great choice if space is at a premium as the arms are slender making the chair narrow and giving it a sleek look. Available in all 4 mechanisms and various backs.



CAMBERLEY

The Camberley is a popular chair because of its simplistic design. Relish those few moments of comfort with a chair that has retro styled arms that is designed to fit in with any suite.

This comfortable, cosy riser recliner is available in the Premier and Deluxe mechanisms with a waterfall back. Its innovative technology brings high quality and great build all for an incredible price. Whilst this chair is only available in a selection of sizes and mechanisms rest assured, we have not compromised on our high quality or expert craftsmanship.

Optional Chair Backs




Waterfall


Standard Fabric Options



Vancouver 21
French Grey fr



Vancouver 21
Mocha fr



Vancouver 21
Praline fr



COMPLETE THE SUITE



3 Seater (Standard) Settee



2 1/2 Seater (Standard) Settee



2 Seater (Standard) Settee



Chair Fixed



Pouffe

Accessories

Arm caps are available for all items and on the rise and recliner chairs they are fastened at the back of the arms as a safety feature. Extra scatter cushions and antimacassars (head throws) are also available. An extra pocket is an available option. A Lithium battery freedom pack with approximately 40-50 manoeuvres before recharging. A foam back spacer is also an option to reduce the seat depth.

Feet Colours (available with fixed furniture)

Natural & Mahogany



Premier Handset



Deluxe Handset



Emergency Battery Back-up (comes as standard)



Safety Cut Off Switch (comes as standard)



Tilt In Space Dual Motor Handset




Tilt In Space Single Motor Handset

*Please contact your local Pride dealer for our full range of accessories and further information


OPTIONAL FABRICS

Other Optional Fabrics


These are non stock items, this will be subject to a manufacturing lead time and increase in price.



Sea Mist




Calico




Ivory




Stone




Caramel




Pebble




Hazelnut




Mist




Willow




Quartz




Wedgwood




Sapphire




Monochrome




Regency Grey



Lead



Pewter



Ash

Shelly Leather Finishes

These are non stock items, this will be subject to a manufacturing lead time and increase in price.



Parchment



Firestone



Tuscany



West



Wood Burner



Crimson



Parlour Blue



Iceblast



Silver Grey



Moon Mist



Silver Birch



Cotton Seed



Cream



Beige



Grove



Almond



Rice Milk



Winter White



White



Seely



Panna



Stone



Angel



Chartreuse



Field Green



Mountain Tree



Pea Green



Thyme Green



Lichen



Horizon



Cherry



Flame Red



Velvet Red



Poppy



Brick Red



Steal



Piping

Shelly Leather Continued

These are non stock items, this will be subject to a manufacturing lead time and increase in price.



Knight



Black



Suffolk Blue



Majolica Blue



Ivory



Somerset
Stone



Basket



Saddle



Spice



Castagna



Deluca



Pebble



Dark Beige



Ash



Owl



Rocking



Jade Green



Forest Green



Golders Green



Sage



Mocca



Dark Chocolate



Anemone



Burgundy



Philly



Blossom



Amethyst



Wineberry



Dark Teal



Bilberry



Deep
Ultramarine



Dark Grape



Havanne



Burnt Oak

Infiniti Faux Leather

These are non stock items, this will be subject to a manufacturing lead time and increase in price.



Caramel



Tan



Espresso



Nugget



Leaf



Racing




Lagoon



Sapphire




Salsa



Granite



Mousse




Anthracite




Bordeaux



Shadow




Taupe



Elephant



Charcoal



Noir



Stone



Camel



Truffle

Unique protection offered

Infiniti conforms to the British Standard for light fastness, tear strength and seam slippage.

This fabric fights against Microbial growth for example: Altinata, Asperigilus Niger, Cladosporium Herbarism, E Coli, Penicillium, Staphylococcus Aureus and M.R.S.A

SPECIFICATIONS

Premier and Deluxe Sizes

	Overall Height	Overall Depth	Overall Width	Seat Height	Seat Depth	Seat Width	Weight Capacity
Premier Grande	1194mm (47")	965mm (38")	915mm (36")	521mm (20 ½")	546mm (21 ½")	533mm (21")	25 st (158.7 kg)
Premier Standard	1143mm (45")	940mm (37")	889mm (35")	495mm (19 ½")	521mm (20 ½")	533mm (21")	25 st (158.7 kg)
Premier Compact	1092mm (43")	915mm (36")	889mm (35")	470mm (18 ½")	495mm (19 ½")	533mm (21")	25 st (158.7 kg)
Premier Petite	1041mm (41")	889mm (35")	889mm (35")	445mm (17 ½")	470mm (18 ½")	533mm (21")	25 st (158.7 kg)
Deluxe Grande	1194mm (47")	965mm (38")	915mm (36")	521mm (20 ½")	546mm (21 ½")	533mm (21")	25 st (158.7 kg)
Deluxe Standard	1143mm (45")	940mm (37")	889mm (35")	495mm (19 ½")	521mm (20 ½")	533mm (21")	25 st (158.7 kg)
Deluxe Compact	1092mm (43")	915mm (36")	889mm (35")	470mm (18 ½")	495mm (19 ½")	533mm (21")	25 st (158.7 kg)
Deluxe Petite	1041mm (41")	889mm (35")	889mm (35")	445mm (17 ½")	470mm (18 ½")	533mm (21")	25 st (158.7 kg)
3 Seater Settee	1117mm (44")	965mm (38")	1828mm (72")	495mm (19 ½")	508mm (20")	1498mm (59")	25 st (158.7 kg)
2 ½ Seater Settee	1117mm (44")	965mm (38")	1575mm (62")	495mm (19 ½")	508mm (20")	1244mm (49")	20 st (127 kg)
2 Seater Settee	1117mm (44")	965mm (38")	1321mm (52")	495mm (19 ½")	508mm (20")	990mm (39")	25 st (158.7 kg)
Chair	1117mm (44")	965mm (38")	810mm (32")	495mm (19 ½")	508mm (20")	482mm (19")	20 st (127 kg)
Pouffe	393mm (15 ½")	533mm (21")	533mm (21")	N/A	N/A	N/A	20 st (127 kg)

Dual and Single Tilt In Space Sizes

	Overall Height	Overall Depth	Overall Width	Seat Height	Seat Depth	Seat Width	Weight Capacity
Dual TIS Grande	1117mm (44")	965mm (38")	915mm (36")	521mm (20 ½")	508mm (20")	559mm (22")	20 st (127 kg)
Dual TIS Standard	1092mm (43")	965mm (38")	864mm (34")	495mm (19 ½")	508mm (20")	508mm (20")	20 st (127 kg)
Dual TIS Petite	1041mm (41")	965mm (38")	864mm (34")	470mm (18 ½")	508mm (20")	508mm (20")	20 st (127 kg)
Single TIS Grande	1117mm (44")	965mm (38")	915mm (36")	521mm (20 ½")	508mm (20")	559mm (22")	20 st (127 kg)
Single TIS Standard	1092mm (43")	965mm (38")	864mm (34")	495mm (19 ½")	508mm (20")	508mm (20")	20 st (127 kg)
Single TIS Petite	1041mm (41")	965mm (38")	864mm (34")	470mm (18 ½")	508mm (20")	508mm (20")	20 st (127 kg)
3 Seater Settee	1092mm (43")	965mm (38")	1854mm (73")	495mm (19 ½")	508mm (20")	1524mm (60")	20 st (127 kg)
2 ½ Seater Settee	1092mm (43")	965mm (38")	1600mm (63")	495mm (19 ½")	508mm (20")	1270mm (50")	20 st (127 kg)
2 Seater Settee	1092mm (43")	965mm (38")	1346mm (53")	495mm (19 ½")	508mm (20")	1016mm (40")	20 st (127 kg)
Chair	1092mm (43")	965mm (38")	864mm (34")	495mm (19 ½")	508mm (20")	508mm (20")	20 st (127 kg)

NEW



WEYMOUTH


Innovative technology meets contemporary style for a recliner that supports you and your busy life. The Weymouth by VivaLift is the perfect addition to your home. The single motor mechanism allows you to fully recline at the touch of a button. Whilst our dual motor mechanism lets you move the back of the chair and footrest independently, helping you find the perfect position. Also letting you fully recline with infinite, lay-flat positioning and experience maximum relaxation.




Mechanism Options:

- Single Motor
- Dual Motor
- Dual Motor with Powered Headrest and Powered Lumbar


Fabric Options:



Merino Fresh



Merino Fawn



Merino Latte

Elevate Your Lifestyle!




Weymouth


Dual Motor lift chair with a powered headrest and powered lumbar support.

The Weymouth dual motor lift chair with a powered headrest and powered lumbar puts you in complete control of your comfort. Not only does this allow you to independently adjust the backrest and footrest, it gives you power to adjust the headrest which is great for getting accurate support to your neck and shoulders, whilst the powered lumbar support offers vital pressure relief to a key part of your back. Helping give you the ultimate relaxation and comfort you deserve. Also with a stylish LED remote with a USB charging point for your phone letting you make the most of your precious time.


Fabric Options:



Merino Fresh



Merino Fawn



Merino Latte




The Wendover

Dual Motor liftchair with luxurious wipe clean faux leather.


Key Features: available in fabric and vinyl finishes, LED control with USB charging point and is a Dual Motor layflat liftchair.

This petite Dual Motor lift chair features two fabric options and a soft touch brown vinyl. Perfect for any living environment The Wendover is a compact electric lift chair, which gently moves you to your perfect position.


Fabric Options:



Merino Fresh



Merino Latte



Brown Vinyl




660 Mini Lounger

A compact chair with limitless comfort and infinite amount of positioning.


Key Features: infinite number of positions, weight capacity of 27 stone and a powerful Dual Motor design.

The 660 Mini Lounger Duet is a Dual Motor lift chair, which features a 27 stone weight capacity and has a waterfall back enabling the upholstery to be easily cleaned. The 660 Mini Lounger is available in three unique fabric options.


Fabric Options:



Merino Fresh



Merino Fawn



Merino Latte




670 Chairbed

Limitless positions available and comfortable waterfall backrest.


Key Features: infinite number of positions, weight capacity of 27 stone, removable easy to clean back rest and a 4-cushion waterfall back.

The 670 Chairbed is a Dual Motor lift chair, which features a 27 stone weight capacity. It has a luxurious waterfall back, enabling the upholstery to be easily cleaned and can be fully folded flat to be transformed into a chairbed. The 670 is available in three unique fabric options.


Fabric Options:



Merino Fresh



Merino Fawn



Merino Latte




LC101

A stylish chair that suits all homes.


Key Features: quiet and smooth 3-position lift system, comfortable T-backrest, full recline and has a powerful Single Motor.

The LC101 comes with a unique T-back design, which makes the LC101 an extremely popular choice within our Single Motor range of lift chairs. Featuring a 3-position lift mechanism, a smooth incline and recline mechanism and available in three luxurious fabric options.


Fabric Options:



Merino Fresh



Merino Fawn



Merino Latte

Specifications

	Overall Width	Seat Depth	Seat Width	Seat to Floor	Top of Back to Seat	Weight Capacity
Weymouth	838mm (33")	508mm (20")	533mm (21")	508mm (20")	635mm (25")	27 st (170 kg)
The Wendover	584mm (23")	444mm (17 ½")	460mm (18")	470mm (18 ½")	584mm (23")	25 st (158 kg)
660	710mm (28")	508mm (20")	508mm (20")	480mm (19")	635mm (25")	27 st (170 kg)
670	762mm (30")	528mm (20 ½")	508mm (20")	508mm (20")	685mm (27")	27 st (170 kg)
LC101	762mm (30")	508mm (20")	508mm (20")	470mm (18 ½")	700mm (27 ½")	23 st (147 kg)

PRIDE[®]

MOBILITY PRODUCTS LTD.

LIVE YOUR BEST[®]

Also from Pride:



Powerchair Range



Large Scooter Range



Travel Scooter Range

Pride Mobility LTD
32 Wedgwood Road
Bicester, Oxon
OX26 4UL, England

Tel: 01869 324600
Email: sales@pride-mobility.co.uk
Web: www.pride-mobility.co.uk



Pride Mobility LTD reserve the right to alter, at any time, without prior notice, the design, specification, packaging or price without the incurring any obligation.