

## Bold, Adventurous and Innovative

The BAJA Wrangler is the new and improves big brother of Prides very own Ranger. The Wrangler is great for rough terrian, with great suspension and an innovative breaking system. Travelling off-road at speeds of up to 8mph has never been easier! Furthermore with optional 100Ah batteries this BAJA powerhouse will keep going all day long at no extra weight.





| Specifications                    |  |
|-----------------------------------|--|
| Scooter                           |  |
| Weight Capacity                   | 32 st 2 lbs (204 kg)                           |
| Overall Length                    | 1556 mm (61.25")                               |
| Overall Width                     | 837 mm (33")                                   |
| Front/Rear Tyres                  | 368 mm (14.5") Solid                           |
| Ground Clearance                  | 140 mm (5.5")                                  |
| Colour(s)                         | Desert Sand, True Blue                         |
| Mobility                          |  |
| Turning Radius                    | 2407 mm (94.75")                               |
| Max Speed <sup>1</sup>            | 8 mph (12.9 kmh)                               |
| Battery                           |  |
| Battery Requirements              | 70 Ah Lead-Acid<br>100 Ah Lead-Acid (optional) |
| Battery Weight                    | 21.75 kg each (70Ah)<br>27.5 kg each (100Ah)   |
| Battery Charger                   | 8-amp Off-Board                                |
| Total Weight<br>Without Batteries | 124.25 kg                                      |

## **Features**



Optional 100 Ah batteries for longer range.



14" Front and rear tyres



Added suspension for a smooth ride even over rough terrain



