



Colours Available:



Rugged, Durable Made For Adventure

Your adventure awaits, and nothing stands in your way. Make movement with the Ranger®, an aggressive, outdoor scooter designed for extreme and rugged conditions. With powerful suspension, bright lighting and dual hydraulic brakes, this beast of a scooter is a powerful force.

Features



Full LED Lights Package



USB Charger Located on the Tiller



LCD Delta Tiller Display



14.5" Chunky Drive Tyres

| Specifications | |
|--|---|
| Scooter | |
| Weight Capacity | 25 st (159 kg) |
| Overall Length ³ | 1556 mm (61.25") |
| Overall Width ³ | 829 mm (32.5") |
| Drive Tyres | 368 mm (14.5") - Pneumatic |
| Ground Clearance ³ | 140 mm (5.5") |
| Colour(s) | Desert Sand, Stone Grey |
| Suspension | Enhanced Front and Rear Suspension |
| Heaviest Piece When Disassembled ³ | 99.25 kg |
| Mobility | |
| Turning Radius ³ | 2407 mm (94.75") |
| Max Speed ¹ | 8 mph (12.8 km/h) |
| Battery | |
| Battery Requirements ³ | 2 x 75 Ah Batteries 2 x 100 Ah Batteries |
| Battery Weight ⁴ | 21.75 kg |
| Battery Charger | 8 - Amp Off-Board |
| Total Weight Without Batteries ^{3, 5} | 124.25 kg |

NOTE:

- 1 Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tyre condition. This specification can be subject to a variance of + or -) 10%.
- 2 Includes standard seat weight.
- 3 Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.
- 4 Battery weight may vary based on battery manufacturer.
- 5 Does not meet airline requirements.

All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.

