

Elevate Your Life To The Next Level!

At the touch of a switch, the Jazzy Air[®] 2 elevates 12 inches in just 11 seconds. With a sleek, modern look, this power chair safely drives at speeds up to 4 mph while elevated, allowing users to enjoy face-to-face social engagement.





Features



Fold Flat Seat



Drives Safely at 4mph While Elevated at 12"



Angle, and Height-Adjustable Headrest

NOTE:

Specifications		NOTE:
Chair		 Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tyre condition. This specification can be subject to a variance of + or -) 10%. Includes standard seat weight. Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.
Weight Capacity	21 st 6 lbs (136 kg)	
Overall Length ³	1067 mm (41.93")	
Overall Width ³	590 mm (23.19")	
Drive Tyres	Drive: 254 mm (10") Flat Free Castors: 152 mm (6") Flat Free	4 Battery weight may vary based on battery manufacturer.
Ground Clearance ³	25 mm (1")	5 Does not meet airline requirements.
Suspension	Active-Trac Suspension With Mid-Wheel 6 Technology	All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of $+$ or -3 %. Dimensions listed are for power
Colour(s)	Citrine yellow, pearly white, ruby red, sapphire blue, black pearl, garnet red, black onyx, pink topaz, silver, tanzanite	
Heaviest Piece When Disassembled ³	82 kg	base only. Overall measurements will vary based on seating and accessory selections.
Power		
Turning Radius ³	528 mm (20.75″)	
Max Speed ¹	4 mph (6 km/h)	
Battery		
Battery Requirements ³	2 x 12v, Deep Cycle U-1 Batteries	
Battery Weight⁴	11.5 kg	
Battery Charger	5 - Amp Off-Board	
Total Weight Without Batteries ^{3, 4}	104.5kg	UK CE



www.pride-mobility.co.uk