



Colours Available:



Style and Comfort At 6 mph

Embrace your potential! The Colt Deluxe 2.0 is fully loaded with style and advanced features, including frosted LED headlights and an integrated storage space with cup holder. Drive at speeds up to 6 mph and see what life has in store for you!

Specifications	
Scooter	
Weight Capacity	28 st 8 lbs (181 kg)
Overall Length	1194 mm (47")
Overall Width	560 mm (22")
Drive Tyres	254 mm Front & Rear (10")- Solid
Ground Clearance	64 mm (2.5")
Heaviest Piece When Disassembled	Front Section 29.5 kg
Colour(s)	Red, Matte Blue, Grey
Suspension	Front & Rear Suspension Front & Rear Suspension
Mobility	
Turning Radius	1334 mm (52.5")
Max Speed	Up to 6.25 mph (10 km/h)
Max. Obstacle Climb	80 mm (3.25")
Drivetrain	24v DC Motor
Dual Breaking System	Regenerative & Electromechanical
Battery	
Battery Requirements	2 x 40 Ah Deep Cycle
Battery Charger	5-A Off-Board
Battery Weight ⁴	14.75 kg Each 40 Ah
Total Weight Without Batteries	68 kg

Features



Comfortable High Back Seat



Integrated USB Charging Point for Device



Pneumatic Tyres With Solid Tyres Available