



Colours Available:



Style and Comfort At 6 mph

Embrace your potential! The Colt Deluxe 2.0 is fully loaded with style and advanced features, including frosted LED headlights and an integrated storage space with cup holder. Drive at speeds up to 6 mph and see what life has in store for you!



Features



Comfortable High Back Seat



Integrated USB Charging Point for Device



Pneumatic Tyres With Solid Tyres Available

Specifications	
Scooter	
Weight Capacity	28 st 8 lbs (181 kg)
Overall Length ³	1200 mm (47.25")
Overall Width ³	563 mm (22.25")
Drive Tyres	269 mm Front & Rear - Solid* or Pneumatic Tyres Available
Ground Clearance ³	64 mm (2.5") at Motor
Heaviest Piece When Disassembled ³	30.25 kg
Colour(s)	Titanium Grey, Ocean Blue, Candy Apple Red
Suspension	Front & Rear Suspension Front & Rear Suspension
Mobility	
Turning Radius ³	1334 mm (52.5")
Max Speed ¹	Up to 6 mph
Max. Obstacle Climb	80 mm (3.25")
Drivetrain	24v DC Motor
Dual Breaking System	Regenerative & Electromechanical
Battery	
Battery Requirements ³	2 x 40 Ah Deep Cycle
Battery Charger	5-Amp Off-Board
Battery Weight ⁴	14.75 kg
Total Weight Without Batteries ^{3, 4}	77.25 kg

NOTE:

- 1 Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tyre condition. This specification can be subject to a variance of + or -) 10%.
- 2 Includes standard seat weight.
- 3 Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.
- 4 Battery weight may vary based on battery manufacturer.
- 5 Does not meet airline requirements.

All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.

