



Colours Available:



Powerful Performance With Exceptional Style

You're active and always on the go. You need a scooter that is up to the task! Appeal to your sporty side with the Colt Executive. Enjoy premium suspension, bright lighting and other sophisticated features, all in one sporty package.

Specifications	
Scooter	
Weight Capacity	28 st 8 lbs (181 kg)
Overall Length	1422 mm (56")
Overall Width	692 mm (27")
Drive Tyres	330 mm Front & Rear (13") - Solid
Ground Clearance	1333 mm (5.25")
Heaviest Piece When Disassembled	Main Frame: 90.3 kg
Colour(s)	Yellow, Black, Grey
Suspension	Front & Rear Suspension
Mobility	
Turning Radius	1905 mm (75.4")
Max Speed	Up to 8 mph (12.5 km/h)
Max. Obstacle Climb	120 mm (4.72")
Drivetrain	24v DC Motor
Dual Breaking System	Regenerative & Electromechanical
Battery	
Battery Requirements	2 x 100 Ah Deep Cycle
Battery Charger	8 A Off-Board
Battery Weight	27.5 kg Each
Total Weight Without Batteries	125.25 kg

Features



Front & Rear LED Lighting Package



Integrated USB Charging Point for Device



Convenient Storage in Tiller