



Colours Available:



# Powerful Performance With Exceptional Style

Enjoy premium suspension, bright lighting, a weight capacity of up to 28st 8lbs, speeds of up to 8mph and other sophisticated features, all in one sporty package.



Specifications	
<b>Scooter</b>	
Weight Capacity	28 st 8 lbs (182 kg)
Overall Length	1422 mm (56")
Overall Width	692 mm (27")
Drive Tyres	Front: 330 mm (13") - Pneumatic Rear: 330 mm (13") - Pneumatic
Heaviest Piece When Disassembled	90.3 kg
Colour(s)	Yellow, Black, Grey
Suspension	Front & Rear Suspension
<b>Mobility</b>	
Turning Radius	1905 mm (75")
Max Speed	8 mph (12.9 km/h)
Max. Obstacle Climb	100 mm
Drivetrain	24v DC Motor
Dual Breaking System	Regenerative & Electromechanical
<b>Battery</b>	
Battery Requirements	2 x 100 Ah Deep Cycle
Battery Charger	8 A Off-Board
Battery Weight	27.5 kg
Total Weight Without Batteries	125 kg

## Features



Front & Rear LED Lighting Package



Convenient Storage in Tiller



REV 02 - 31/01/2023

