



Colours Available:



Powerful Performance With Exceptional Style

You're active and always on the go. You need a scooter that is up to the task! Appeal to your sporty side with the Colt Executive. Enjoy premium suspension, bright lighting and other sophisticated features, all in one sporty package.

Features



Front & Rear LED Lighting Package



Integrated USB Charging Point for Device



Convenient Storage in Tiller

Specifications	
Scooter	
Weight Capacity	28 st 8 lbs (182 kg)
Overall Length ³	1428 mm (56.25")
Overall Width ³	686 mm (27")
Drive Tyres	330 mm Front & Rear - Pneumatic
Ground Clearance ³	1333 mm (5.25")
Heaviest Piece When Disassembled ³	Main Frame: 90.3 kg
Colour(s)	Titanium Grey, Black & Yellow
Suspension	Front & Rear Suspension
Mobility	
Turning Radius ³	1973 mm (76.25")
Max Speed ¹	Up to 8 mph
Max. Obstacle Climb	120 mm (4.72")
Drivetrain	24v DC Motor
Dual Breaking System	Regenerative & Electromechanical
Battery	
Battery Requirements ³	2 x 100 Ah Deep Cycle
Battery Charger	8-Amp Off-Board
Battery Weight ⁴	27.4 kg
Total Weight Without Batteries ^{3,4}	117 kg

NOTE:

- 1 Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tyre condition. This specification can be subject to a variance of + or -) 10%.
- 2 Includes standard seat weight.
- 3 Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.
- 4 Battery weight may vary based on battery manufacturer.
- 5 Does not meet airline requirements.

All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.

