



Colours Available:



Versatile and Stylish Know No Bounds

Enjoy the simple and disassembly Colt Plus. The bright LED headlight package and convenient front storage basket means anything is possible with the Colt Plus!

Specifications	
Scooter	
Weight Capacity	25 st (159 kg)
Overall Length	1194 mm (47")
Overall Width	560 mm (22")
Drive Tyres	254 mm Front & Rear Pneumatic
Ground Clearance	57 mm (2.25")
Heaviest Piece When Disassembled	29.75 kg
Colour(s)	Red
Suspension	Limited
Mobility	
Turning Radius	1370 mm (54")
Max Speed	Up to 4 mph (6 km/h)
Max. Obstacle Climb	78.7 mm (3.1")
Drivetrain	24v DC Motor
Dual Breaking System	Regenerative & Electromechanical
Battery	
Battery Requirements	2 x 32 Ah Deep - Cycle
Battery Charger	4-Amp Off-Board
Total Weight Without Batteries	64.5 kg

Features



Convenient Storage in Tiller



Exclusive & Stylish Low Profile Tyre



High Visibility Front and Rear LED Lighting