



Colours Available:



Versatile and Stylish Know No Bounds

Enjoy the simple and disassembly Colt Plus that comes with a bright LED headlight package as well as a convenient front storage basket.

Specifications	
Scooter	
Weight Capacity	25 st (159 kg)
Overall Length	1194 mm (47")
Overall Width	560 mm (22")
Drive Tyres	Front: 254 mm (10") Pneumatic Rear: 254 mm (10") Pneumatic
Ground Clearance	57 mm (2.25")
Heaviest Piece When Disassembled	29.75 kg
Colour(s)	Red
Suspension	Limited
Mobility	
Turning Radius	1370 mm (54")
Max Speed	4 mph (6.4 km/h)
Max. Obstacle Climb	50 mm
Drivetrain	24v DC Motor
Dual Breaking System	Regenerative & Electromechanical
Battery	
Battery Requirements	2 x 30 Ah (standard)
Battery Charger	4-Amp Off-Board
Total Weight Without Batteries	64.5 kg

Features



Convenient Storage in Tiller



Exclusive & Stylish Low Profile Tyre



High Visibility Front and Rear LED Lighting