



Colours Available:



Candy Apple Red

Versatile and Stylish Know No Bounds

Enjoy the simple and disassembly Colt Plus. The bright LED headlight package and convenient front storage basket means anything is possible with the Colt Plus!

Features



Convenient Storage in Tiller



Exclusive & Stylish Low Profile Tyre



High Visibility Front and Rear LED Lighting

Specifications	
Scooter	
Weight Capacity	25 st (158.7 kg)
Overall Length ³	1190 mm (47")
Overall Width ³	560 mm (22")
Drive Tyres	254 mm Front & Rear - Pneumatic
Ground Clearance ³	57 mm (2.25")
Heaviest Piece When Disassembled ³	29.7 kg
Colour(s)	Candy Apple Red
Suspension	Limited
Mobility	
Turning Radius ³	1370 mm (54")
Max Speed ¹	Up to 4 mph
Max. Obstacle Climb	78.7 mm (3.1")
Drivetrain	24v DC Motor
Dual Breaking System	Regenerative & Electromechanical
Battery	
Battery Requirements ³	2 x 32 Ah Deep Cycle
Battery Charger	4-Amp Off-Board
Total Weight Without Batteries ^{3, 4}	62.5 kg

NOTE:

- 1 Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tyre condition. This specification can be subject to a variance of + or -) 10%.
- 2 Includes standard seat weight.
- 3 Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.
- 4 Battery weight may vary based on battery manufacturer.
- 5 Does not meet airline requirements.

All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.

