

Specifications	
Scooter	
Weight Capacity	25 st (159 kg)
Overall Length	1359 mm (53")
Overall Width	635 mm (25")
Overall Height	1346 mm (53")
Drive Tyres	Front & Rear - 330 mm (13")
Ground Clearance	94 mm (3.7")
Heaviest Piece When Disassembled	20.11 kg
Colour(s)	Silver & Black
Suspension	Front & Rear Suspension
Standard Seating	Type: CRS Deluxe Weight: 19 kg (43 lbs) Material: Black Vinyl Width: 457 mm (18") Depth: 457 mm (18")
Mobility	
Turning Radius	1606 mm (63.13")
Max Speed	Up to 8 mph(12.87km/h)
Max. Obstacle Climb	100 mm (4")
Drivetrain	Rear - Wheel Drive Sealed Transaxle 24v DC Motor
Dual Breaking System	Regenerative & Electromechanical
Maximum Obstacle Climbing Ability	100 mm (4")
Battery	
Battery Requirements	75 Ah 2 x 12v Deep - Cycle
Battery Weight	22.11kg Each
Battery Charger	8 A Off-Board Charger
Total Weight Without Batteries	96 kg

Features



High Visibility Front and Rear LED Lighting



Indicators



13" Alloy Style Pneumatic Tyres