



Colours Available:



Robust Features With Supreme Quality


The Colt Pursuit ES13 drives at speeds of up to 8 mph. Equipped with a regenerative braking system, a comfortable high-back seat and complete LED lighting, the ES13 meets your every need.

Specifications	
Scooter	
Weight Capacity	25 st (159 kg)
Overall Length	1355 mm (53.5")
Overall Width	635 mm (25")
Drive Tyres	Front: 330 mm (13") pneumatic Rear: 330 mm (13") pneumatic
Ground Clearance	94 mm (3.7")
Heaviest Piece When Disassembled	20.11 kg
Colour(s)	Silver & Black
Suspension	Front & Rear Suspension
Standard Seating	Type: CRS Deluxe Weight: 19 kg (43 lbs) Material: Black Vinyl Width: 457 mm (18") Depth: 457 mm (18")
Mobility	
Turning Radius	1606 mm (63.13")
Max Speed	8 mph (12.9 km/h)
Max. Obstacle Climb	100 mm
Drivetrain	Rear - Wheel Drive Sealed Transaxle 24v DC Motor
Dual Breaking System	Regenerative & Electromechanical
Battery	
Battery Requirements	2 x 50 Ah
Battery Weight	22.11kg
Battery Charger	8 A Off-Board Charger
Total Weight Without Batteries	96 kg

Features



High Visibility Front and Rear LED Lighting



Indicators



13" Alloy Style Pneumatic Tyres



REV 02 - 31/01/2023

