



Colours Available:



# Robust Features With Supreme Quality

Do you have the need for speed? The Colt Pursuit ES13 drives at speeds up to 8 mph. Equipped with a regenerative braking system, a comfortable high-back seat and complete LED lighting, the ES13 appeals to the adventurer within you!

## Features



High Visibility Front and Rear LED Lighting



Indicators



13" Alloy Style Pneumatic Tyres

Specifications	
<b>Scooter</b>	
Weight Capacity	25 st (159 kg)
Overall Length <sup>3</sup>	1359 mm (53.5")
Overall Width <sup>3</sup>	635 mm (25")
Drive Tyres	330 mm Front & Rear - Pneumatic
Ground Clearance <sup>3</sup>	89 mm (3.5")
Heaviest Piece When Disassembled <sup>3</sup>	19 kg
Colour(s)	Silver & Black
Suspension	Front & Rear Suspension
<b>Mobility</b>	
Turning Radius <sup>3</sup>	1740 mm (68.5")
Max Speed <sup>1</sup>	Up to 8 mph
Max. Obstacle Climb	100 mm (4")
Drivetrain	24v DC Motor
Dual Breaking System	Regenerative & Electromechanical
<b>Battery</b>	
Battery Requirements <sup>3</sup>	2 x 75 Ah Deep Cycle
Battery Weight <sup>4</sup>	24 kg
Battery Charger	8-amp Off-Board
Total Weight Without Batteries <sup>3, 4</sup>	73 kg

### NOTE:

- 1 Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tyre condition. This specification can be subject to a variance of + or - 10%.
- 2 Includes standard seat weight.
- 3 Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.
- 4 Battery weight may vary based on battery manufacturer.
- 5 Does not meet airline requirements.

All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3 %. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.