



Colours Available:



Sporty Style With a Advanced Suspension System

Embrace your potential! The Colt Sport is fully loaded with sporty style and advanced features, including frosted LED headlights and an advanced suspension system. Drive at speeds up to 8 mph and see what life has in store for you!

Specifications	
Scooter	
Weight Capacity	25 st (159 kg)
Overall Length	1194 mm (47")
Overall Width	657 mm (25.9")
Drive Tyres	254 mm Front & Rear - (10") Solid
Ground Clearance	57 mm (2.25")
Heaviest Piece When Disassembled	29.5 kg
Colour(s)	Blue & Red
Suspension	Front & Rear Suspension
Mobility	
Turning Radius	1334 mm (52.5")
Max Speed	Up to 8 mph (12.9 km/h)
Max. Obstacle Climb	50 mm (2")
Drivetrain	24v DC Motor
Dual Breaking System	Regenerative & Electromechanical
Battery	
Battery Requirements	2 x 40 Ah Deep Cycle.
Battery Weight	11 kg
Battery Charger	Off-Board
Total Weight Without Batteries	68 kg

Features



Indicators



13" Alloy Style Pneumatic Tyres



Comfortable High Back Seat
With Armrests