



Colours Available:



# Sporty Style With a Advanced Suspension System

Embrace your potential! The Colt Sport is fully loaded with sporty style and advanced features, including frosted LED headlights and an advanced suspension system. Drive at speeds up to 8 mph and see what life has in store for you!



Features



Comfortable High Back Seat With Armrests



Indicators



13" Alloy Style Pneumatic Tyres

Specifications	
<b>Scooter</b>	
Weight Capacity	25 st (159 kg)
Overall Length <sup>3</sup>	1194 mm (47")
Overall Width <sup>3</sup>	657 mm (25.9")
Drive Tyres	254 mm (10") Front & Rear - Solid
Ground Clearance <sup>3</sup>	57 mm (2.25")
Heaviest Piece When Disassembled <sup>3</sup>	29.5 kg
Colour(s)	Red & Blue
Suspension	Front & Rear Suspension
<b>Mobility</b>	
Turning Radius <sup>3</sup>	1333.5 mm (52.5")
Max Speed <sup>1</sup>	Up to 8 mph
Max. Obstacle Climb	50mm (2")
Drivetrain	24v DC Motor
Dual Breaking System	Regenerative & Electromechanical
<b>Battery</b>	
Battery Requirements <sup>3</sup>	2 x 40 Ah Deep Cycle
Battery Weight <sup>4</sup>	14.75 kg
Battery Charger	Off-Board
Total Weight Without Batteries <sup>3, 4</sup>	68 kg

NOTE:

- 1 Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tyre condition. This specification can be subject to a variance of + or -) 10%.
- 2 Includes standard seat weight.
- 3 Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.
- 4 Battery weight may vary based on battery manufacturer.
- 5 Does not meet airline requirements.

All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.

