

Fusion



Colours Available:



Innovative, Compact Design

Highly Manoeuvrable In Tight Spaces

The Fusion rear-wheel drive power chair delivers aggressive, terrain-gripping torque with a maximum speed of up to 4 mph (6.43 km/h).

The innovative, compact design of the Fusion makes it highly manoeuvrable in tight spaces without sacrificing traditional rear-wheel drive performance.

Specifications	
Chair	
Weight Capacity	21 st 6 lbs (136 kg)
Overall Length	890 mm (35")
Overall Width	635 mm (25")
Drive Tyres	355 mm (14") Pneumatic
Ground Clearance	110 mm (4.5")
Colour(s)	Silver
Suspension	Rear
Castors	230 mm (9") Solid Front
Power	
Turning Radius	530 mm (21")
Max Speed	Up to 6 mph (10 km/h)
Drivetrain	Two-Motor, Rear-Wheel Drive
Dual Breaking System	Intelligent Braking (Regenerative, Electromagnetic Disc Park Brake)
Battery	
Battery Requirements	2 x 50 Ah Deep Cycle
Battery Charger	5A Off-Board
Battery weight	17 kg
Total Weight Without Batteries	74 kg

Features



Swing-Away Joystick



Handy Footstool



Individual Freewheel Levers