

Fusion



Colours Available:



Innovative, Compact Design

Highly Manoeuvrable In Tight Spaces

The Fusion rear-wheel drive power chair delivers aggressive, terrain-gripping torque with a maximum speed of up to 4 mph (6.4 km/h).

The innovative, compact design of the Fusion makes it highly manoeuvrable in tight spaces without sacrificing traditional rear-wheel drive performance.

Specifications	
Chair	
Weight Capacity	21 st 6 lbs (136 kg)
Overall Length	890 mm (35")
Overall Width	635 mm (25")
Front/Rear Tyres	Front: 230mm (9") Pneumatic Rear: 355 mm (21") Pneumatic
Ground Clearance	110 mm (4.5")
Colour(s)	Silver
Suspension	Rear
Power	
Turning Radius	530 mm (21")
Max Speed	5 mph (8 km/h)
Drivetrain	Two-Motor, Rear-Wheel Drive
Dual Breaking System	Intelligent Braking (Regenerative, Electromagnetic Disc Park Brake)
Battery	
Battery Requirements	2 x 50 Ah Deep Cycle
Battery Charger	5A Off-Board
Battery weight	17 kg
Total Weight Without Batteries	74 kg

Features



Swing-Away Joystick



Handy Footstool



Individual Freewheel Levers



REV 02 - 30/01/2023

