

Jazzy[®]
600 ES



Colours Available:




Extraordinary Manoeuvrability With a Flawless Design

The Jazzy[®] 600 ES has 14-inch drive wheels and independent front casters that handle terrain changes easily. Enjoy aggressive precision, a tight turning radius and a flawless design.


Features



Fold Flat Seat



Extra Storage



Angle, Height and Depth-Adjustable Armrest

Specifications	
Chair	
Weight Capacity	21 st 6 lbs (136 kg)
Overall Length ³	900.2mm
Overall Width ³	622 mm (24.5")
Drive Tyres	Drive: 3556 mm (14") Solid Castors: 1524 mm (6") Solid
Ground Clearance ³	70 mm (2.75") at Battery Tray
Suspension	ATX Suspension (Active-Trac With Extra Stability)
Colour(s)	Black, Blue & Red
Power	
Turning Radius ³	521 mm (20.5") With Front Riggings
Max Speed ¹	4 mph (6 km/h)
Max. Obstacle Climb	70 mm (2.75")
Battery	
Battery Requirements ³	2 x 12v, Deep Cycle NF-22 Batteries
Battery Weight ⁴	16.9 kg each
Battery Charger	8 - Amp Off-Board
Total Weight Without Batteries ^{3, 4}	56.1 kg

NOTE:

- 1 Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tyre condition. This specification can be subject to a variance of + or -) 10%.
- 2 Includes standard seat weight.
- 3 Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.
- 4 Battery weight may vary based on battery manufacturer.
- 5 Does not meet airline requirements.

All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3 %. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.

