

Jazzy[®] SELECT



Colours Available:



Highly Manoeuvrable Loaded With Performance And Convenience

The highly manoeuvrable Jazzy[®] Select delivers the advantages of in-line motor technology for enhanced efficiency, torque, range, and performance. Plus, it's loaded with performance and convenience features, like standard Active-Trac[®] Suspension, that make it easy to enjoy.

Specifications	
Chair	
Weight Capacity	21 st 6 lbs (136 kg)
Max Speed	4 mph
Overall Length	908 mm (35.75")
Overall Width	603 mm (23.75")
Drive Tyres	254 mm (10")
Ground Clearance	47.6 mm (1.875") at Midframe
Suspension	Active-Trac®
Colour(s)	Blue & Red
Power	
Turning Radius	508 mm (35.75")
Drivetrain	Two - Motor, Mid - Wheel Drive
Dual Breaking System	"Intelligent Braking" Electronic Regenerative, Disc Park Brake
Battery	
Battery Requirements	2 x 30 Ah, Deep Cycle
Battery Weight	11.11 kg
Battery Charger	5 - Amp Off-Board

Features



Flip Out Foot Board



Intuitive Controller is Easy to Operate



Active-Trac® Suspension