

# Jazzy<sup>®</sup> SELECT



Colours Available:



## Highly Manoeuvrable Loaded With Performance And Convenience

The highly manoeuvrable Jazzy<sup>®</sup> Select delivers the advantages of in-line motor technology for enhanced efficiency, torque, range, and performance. Plus, it's loaded with performance and convenience features, like standard Active-Trac<sup>®</sup> Suspension, that make it easy to enjoy.

Specifications	
<b>Chair</b>	
Weight Capacity	21 st 6 lbs (136 kg)
Max Speed	4 mph (6.4km/h)
Overall Length	1118 mm (44")
Overall Width	603 mm (23.75")
Front/Rear Tyres	Front: 146mm (5.75") solid Rear: 254 mm (10") solid
Ground Clearance	19 mm (0.75")
Suspension	Active-Trac®
Colour(s)	Blue & Red
<b>Power</b>	
Turning Radius	571 mm (22.5")
Drivetrain	Two - Motor, Mid - Wheel Drive
Dual Breaking System	"Intelligent Braking" Electronic Regenerative, Disc Park Brake
<b>Battery</b>	
Battery Requirements	2 x 30 Ah
Battery Weight	11.5 kg
Battery Charger	3.5 - Amp Off-Board

## Features



Flip Out Foot Board



Intuitive Controller is Easy to Operate



Active-Trac® Suspension



REV 02 - 30/01/2023

